

# Empirical Formula of $Cu_xCl_y$

NEXT: LAB  $\rightarrow$  Reactivity of Metals

## (a) Empirical Formula Determination (Mole Ratios)



$$Mg = 0.353g$$

$$Mg_x O_y = 0.585g$$

$$0.353g Mg \times \frac{1 \text{ mol } Mg}{24.3g Mg} = 1.45 \times 10^{-2} \text{ mol } Mg$$

$$(0.585 - 0.353)g O \times \frac{1 \text{ mol } O}{16g O} = 1.45 \times 10^{-2} \text{ mol } O$$

$$\frac{1.45 \times 10^{-2} \text{ mol } Mg}{1.45 \times 10^{-2} \text{ mol } Mg} = 1 \quad \frac{1.45 \times 10^{-2} \text{ mol } Mg}{1.45 \times 10^{-2} \text{ mol } O} = 1 \leftarrow \text{smallest value.}$$

$$x = 1 \quad y = 1$$

$$\text{Percent yield } Cu = \frac{\text{Actual Yield } Cu}{\text{Theo. Yield } Cu} \times 100$$

$$\text{Theo. yield } Cu = \boxed{\phantom{00}}g Zn \times \frac{1 \text{ mol } Zn}{65.39g Zn} \cdot \frac{1 \text{ mol } Cu}{1 \text{ mol } Zn} \cdot \frac{63.55g Cu}{1 \text{ mol } Cu} = \boxed{\phantom{00}}g Cu$$

- slip Zn in  $Cu_xCl_y$  sh.
- scrape off Cu from Zn with glass rod, carefully with no spills.
- Take plastic coating off.
- Setup filter + flask + buchner funnel, + aspirator.
- $Cu_xCl_y$  [ ] =
- Reaction flask over white paper.
- Discard Cu in container, Zn sh in sinks, solid Zn in container